

IZBJELJIVANJE Er:YAG LASEROM

IZV.prof.dr.sc. Ivona Bago

Energija: 40 mJ

Frekvencija: 15 Hz

VLP

3 x 20 s svaki zub



Laserski "handpiece" za izbjeljivanje



Er:YAG (Fotona, LightWalker AT)
Opalescence Boost (40% vodíkov peroksid)



Opalescence Boost (40% vodíkov peroksid)



LASERIRANJE

MOD: VLP

40-60 mJ

10-30 Hz

3 x 20 s po zubu



PARAMETRI LASERA

MOD: VLP

40-60 mJ

10-30 Hz

3 x 20 s po zubu









